



The Cardiac Troponin T Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of human cardiac Troponin T in whole blood, serum or plasma as an aid in the diagnosis of myocardial infarction (MI).

Product Features

- Fast results
- Easy visual interpretation
- Simple operation, no equipment required
- High accuracy

Specimen	WB/S/P
Certificate	CE
Reading Time	10 minutes
Pack	10 T
Storage Temperature	2-30°C
Shelf Life	2 Years
Sensitivity	98.6%
Specificity	99.5%
Accuracy	99.4%
Cut-Off	0.5 ng/mL

PRODUCT SPECIFICATION

Ordering Information

Cat. No.	CTNT-402
Product	Cardiac Troponin T (cTnT) Rapid Test
Specimen	WB/S/P
Pack	10 T
Cut-Off	0.5 ng/mL

Tel: 0086-571-56267891

E-mail: info@alltests.com.cn

Add: 550 Yin Hai street, Hangzhou Economic and Technologic Development Area 310018, P.R China